

Onze infuuspompen



CADD Solis VIP



CADD Legacy ONE



Bodyguard 323 Color Vision



Medcaptain MP60



Medcaptain HP60



Medcaptain MP30



Medcaptain HP30



Crono S-pid 30



Crono S-pid 50



Crono S-pid 100



Elastomeer

| Gebruik | Mobiel | Mobiel | Mobiel | Stationair | Stationair | Stationair | Stationair | Mobiel | Mobiel | Mobiel | Mobiel |
|---------|--------|--------|--------|------------|------------|------------|------------|--------|--------|--------|--------|
|---------|--------|--------|--------|------------|------------|------------|------------|--------|--------|--------|--------|

Toedieningswijze

| | | | | | | | | | | | |
|-----------|---|---|---|---|---|---|---|---|---|---|---|
| Subcutaan | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ |
| Perifeer | √ | √ | √ | √ | √ | √ | √ | | | | √ |
| CVC | √ | √ | √ | √ | √ | √ | √ | | | | √ |
| PAC | √ | √ | √ | √ | √ | √ | √ | | | | √ |
| PICC | √ | √ | √ | √ | √ | √ | √ | | | | √ |
| Epiduraal | √ | | | | √ | | √ | | | | |
| Spinaal | √ | | | | √ | | √ | | | | |
| Midline | √ | √ | √ | √ | √ | √ | √ | | | | √ |

Doseringsfrequentie

| | | | | | | | | | | | |
|--------------------|---------|---------|---------|----------|----------|----------|----------|---------|---------|---------|---------|
| Bolus | √ | | √ | √ | √ | √ | √ | √ | √ | √ | |
| Continu ml/uur | √ | | √ | √ | √ | √ | √ | √ | √ | √ | |
| Continu ml/ 24 uur | | √ | | √ | √ | | √ | | | | |
| Intermitterend | √ | | √ | √ | √ | √ | √ | √ | √ | √ | |
| PCA | √ | | √* | √ | √ | √ | √ | | | | |
| Stappenprogramma | √ | | √ | √ | √ | √ | √ | √ | √ | √ | |
| Min. Loopsnelheid | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,1ml/u | 0,2ml/u |
| Max. Loopsnelheid | 50ml/u | 125ml/u | 500ml/u | 1200ml/u | 2000ml/u | 2000ml/u | 2200ml/u | 75ml/u | 100ml/u | 200ml/u | |

Therapieën

| | | | | | | | | | | | |
|--|---|---|---|---|---|---|---|---|---|---|---|
| Medicatie Continu | √ | √ | √ | | | | | | | | √ |
| Medicatie Gift | | | √ | √ | √ | √ | √ | | | | √ |
| Immunoglobulinen Intraveneus | | | √ | | √ | | | | | | |
| Immunoglobulinen Subcutaan | | | | | | | | √ | √ | √ | |
| Gefaciliteerd Immunoglobulinen Subcutaan | | | √ | | | | | | | | |
| Totale Parentale Voeding | √ | | √ | | | | | | | | |
| Pijn en Sedatie | √ | | √ | | | | √ | √ | | | |
| Bloedtransfusie | | | | | √ | | | | | | |
| Oncologie Continu | √ | | | | | | | | | | |
| Oncologie Gift | | | √ | √ | √ | | | | | | |

Alarmfuncties

| | | | | | | | | | | | |
|------------------------------|---|---|---|---|---|---|---|---|---|---|---|
| Batterij/ accu leeg | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ |
| Infusie gereed/ pomp gestopt | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Lucht in de lijn | √ | √ | √ | √ | √ | | | | | | |
| Oclusiedruk | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ |
| Reservoir leeg | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ |
| Systeemfout | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ |

Doseringshoeveelheid

| | | | | | | | | | | | |
|-----------------|---|---|---|---|---|--|---|---|---|---|------------------|
| 50 ml Cassette | √ | √ | √ | | | | | | | | |
| 100 ml Cassette | √ | √ | √ | | | | | | | | |
| 250 ml Cassette | √ | √ | √ | | | | | | | | |
| Infuuszak | √ | √ | √ | √ | √ | | | | | | |
| 5 ml Spuit | | | | | | | | √ | | | |
| 10 ml Spuit | | | | | | | √ | √ | | | |
| 20 ml Spuit | | | | | | | √ | √ | | | |
| 30 ml Spuit | | | | | | | √ | √ | √ | | |
| 50/60 ml Spuit | | | | | | | √ | √ | | √ | |
| Fles | √ | √ | √ | √ | √ | | | | | | |
| Elastomeer | | | | | | | | | | | 50, 100 en 250ml |